

# Calendar Home Page (Medical Cart)

1

Search for and join meetings with your contacts or rooms, or dial out to H.323/SIP endpoints

2

Adjust speaker volume, view enabled devices, update network settings, and view account info

3

View a custom Help page with contact information and important links for troubleshooting

4

Click to enter a meeting ID from a previously created calendar appointment

5

View the date, time, and name of your room


The screenshot displays the VidyoConnect Medical Cart interface. On the left is a vertical navigation menu with icons and labels: Home (house icon), Search (magnifying glass icon), Contacts (two people icon), Rooms (square icon with a dot), Dial out (grid icon), Settings (gear icon), and Help (question mark icon). Yellow callout circles with numbers 1 through 5 are placed over the interface. Callout 1 points to the Search icon. Callout 2 points to the Settings icon. Callout 3 points to the Help icon. Callout 4 points to a green 'Join by Meeting ID' button within a dark overlay that says 'Please enter meeting ID to join a conference'. Callout 5 points to a status card in the top right corner showing 'Tuesday, November 28 6:40 PM' and 'Valley Hospital Connected' with a profile picture. The main content area shows a 'Meet with Patient' card with 'Next call: 3 PM to 3:20 PM', 'Invitees' and 'Join' buttons, and a partially visible 'Call: 5 PM to 5:20 PM' card below it.

# Modify Settings

1 Click to view Settings

2 Click to select Settings options

3 Click the down arrow to scroll down for more Settings (Exit and Reboot)



# Meet Now

- 1 Subject of the meeting
- 2 Meeting start and end times
- 3 Click to view meeting invitees
- 4 Click to join the meeting
- 5 View next call

The screenshot shows the VidyoConnect mobile application interface. On the left is a vertical navigation menu with icons and labels: Home (house icon), Search (magnifying glass icon), Contacts (two people icon), Rooms (square with dot icon), Dial out (grid of dots icon), Settings (gear icon), and Help (question mark icon). The main content area features a dark overlay at the top with the text "Please enter meeting ID to join a conference" and a green "Join by Meeting ID" button. Below this is a white card titled "Meet with Patient" which displays "Next call: 3 PM to 3:20 PM" and two buttons: "Invitees" (blue) and "Join" (green). At the bottom, another white card shows "Call: 5 PM to 5:20 PM". A top right status bar displays "Tuesday, November 28 6:40 PM" and a profile for "Valley Hospital Connected" with a circular photo of two people. Five yellow circular callouts with numbers 1 through 5 are connected by dotted lines to specific elements: 1 points to the "Meet with Patient" title, 2 points to the "Next call" time, 3 points to the "Invitees" button, 4 points to the "Join" button, and 5 points to the "Call: 5 PM to 5:20 PM" text.

# Join a Call by Meeting ID

1

Copy the Meeting ID from a Google or Exchange calendar appointment

2

Click the Join by Meeting ID button and then enter/paste the Meeting ID

3

Join the meeting

The screenshot displays the VidyConnect mobile application interface. On the left is a vertical navigation menu with icons and labels for Home, Search, Contacts, Rooms, Dial out, Settings, and Help. The main content area shows a 'Meet with Patient' card with the text 'Next call: 3 PM to 3:20 PM'. A dark overlay box is positioned over the 'Join by Meeting ID' button, which is highlighted with a yellow circle and the number '2'. Above this button, the text reads 'Please enter meeting ID to join a conference'. In the top right corner, a status bar shows the date 'Tuesday, November 28' and the time '6:40 PM', along with a profile picture and the text 'Valley Hospital Connected'. Below the main card, another card is partially visible with the text 'Call: 5 PM to 5:20 PM'.

# Search to Add People, Rooms, or H.323/SIP Endpoints

- 1 Click to open the Search page
- 2 Use the VidyoControl app or the IR Remote with the keypad to input numbers and text
- 3 View alphabetized results of people, rooms, or endpoints
- 4 Use the arrow keys to navigate and scroll left and right
- 5 Click to select the contact, room, or endpoint

The screenshot shows the VidyoConnect search interface. On the left is a sidebar with navigation icons: Home, Search, Contacts, Rooms, Dial out, Settings, and Help. The main area features a search input field containing 'Dr', a numeric keypad, and a list of search results. Callout 1 points to the Search icon in the sidebar. Callout 2 points to the search input field. Callout 3 points to the search results. Callout 4 points to the horizontal scroll bar above the results. Callout 5 points to the Dr. Collum result.

a	b	c	d	e	f
g	h	i	j	k	l
m	n	<b>o</b>	p	q	r
s	t	u	v	w	x
y	z	0	1	2	3
4	5	6	7	8	9
@	*	#	-	-	/
,	.	:	⏪	ABC	⏩
Clear all					

Search results:

- Dr. Anderson (Person)
- Dr. Anderson's Room (Room)
- Dr. Collum's Room (Room)
- Dr. Brown (Person)
- Dr. Collum (Person)
- Dr. Davidson (Person)


# Add Favorite People, Rooms, or H.323/SIP Endpoints

- 1 Click to add this contact to your favorites
- 2 Click to turn Mic on or off
- 3 Click to turn Camera on or off
- 4 Click to call the contact
- 5 Click to return to previous page

The screenshot shows the VidyoConnect interface with a sidebar on the left and a main content area. The sidebar contains the following items from top to bottom: VidyoConnect logo, Home (house icon), Search (magnifying glass icon), Contacts (two people icon), Rooms (door icon), Dial out (grid icon), Settings (gear icon), and Help (question mark icon). The main content area displays a contact card for 'Dr. Anderson's Room' with a circular profile picture of a stethoscope and a heart. Below the name are several controls: a green 'Join' button, a 'Mic on' button with a microphone icon, a 'Camera on' button with a camera icon, a grey 'Back' button, and a blue 'Add to Favorites' button. Numbered callouts (1-5) are placed around the interface to indicate the steps: 1 points to 'Add to Favorites', 2 points to 'Mic on', 3 points to 'Camera on', 4 points to 'Join', and 5 points to 'Back'.

# Contacts Page with Search

- 1 Click to open Contacts page
- 2 Click to search for contacts
- 3 Use your arrow keys to navigate and scroll left and right



# Search for Contacts

1

Use the VidyoControl app or the IR Remote with the keypad to input numbers and text to enter a contact name

2

View alphabetized results of people (contacts)

3

Use the arrow keys to navigate and scroll left and right

4

Click to select the contact to join the meeting

The screenshot shows the VidyoConnect app interface. On the left is a sidebar with icons for Home, Search, Contacts, Rooms, Dial out, Settings, and Help. The main area is titled 'Frequent contacts (107)' and 'Favorites (61)'. A search bar at the top contains the text 'And'. Below it is a keypad with letters a-z, numbers 0-9, and symbols like @, \*, #, -, and /. The letter 'o' is highlighted in blue. To the right of the keypad, two contact cards are visible: 'Dr. Anderson' and 'Dr. Ander'. A yellow circle with the number '1' points to the search bar. A yellow circle with the number '2' points to the 'Dr. Anderson' contact card. A yellow circle with the number '3' is at the top right with a double-headed arrow. A yellow circle with the number '4' points to the 'Dr. Ander' contact card.

VidyoConnect™

Home

Search

Contacts

Rooms

Dial out

Settings

Help

Frequent contacts (107)

Favorites (61)

And

a	b	c	d	e	f
g	h	i	j	k	l
m	n	o	p	q	r
s	t	u	v	w	x
y	z	0	1	2	3
4	5	6	7	8	9
@	*	#	-	-	/
,	.	:	⏪	ABC	⏩

Clear all


Dr. Anderson

Dr. Ander




# Join Meetings with Contacts

- 1 Click to turn Mic on or off
- 2 Click to turn Camera on or off
- 3 Click to call the contact (button will be grayed out if offline or busy)
- 4 Click to remove a favorite contact
- 5 Click to return to the previous page a favorite contact



# Rooms Page with Search

- 1 Click to open Search page
- 2 Click to search for Rooms
- 3 Use the arrow keys to navigate and scroll left and right



# Search for Rooms

1

Use the VidyoControl app or the IR Remote with the keypad to input numbers and text

2

View alphabetized sorted results of rooms to join a meeting

3

Use the arrow keys to navigate and scroll left and right

4

Click to select the room to join a meeting

VidyoConnect™

Home

Search

Contacts

Rooms

Dial out

Settings

Help

Frequent rooms (107)

Favorites (61)

Office

a	b	c	d	e	f
g	h	i	j	k	l
m	n	<b>o</b>	p	q	r
s	t	u	v	w	x
y	z	0	1	2	3
4	5	6	7	8	9
@	*	#	-	-	/
,	.	:	⏪	ABC	⏩

Clear all

Dr. Anderson's Office

Dr. Collum's Office

Dr. Davidson's Office

Dr. Office #3

Dr. Brown's Office

Dr. Office #1

Dr. Office #2

Dr. Office #4

# Join Meetings in Rooms

- 1 Click to turn Mic on or off
- 2 Click to turn Camera on or off
- 3 Click to join the meeting
- 4 Click to remove a favorite room
- 5 Click to return to previous page

The screenshot displays the VidyoConnect interface for a meeting titled "Lab Room". On the left is a vertical navigation menu with icons and labels: VidyoConnect™, Home, Search, Contacts, Rooms (highlighted in green), Dial out, Settings, and Help. The main content area shows a circular profile picture of a person, the room name "Lab Room", and a green "Join" button. Below the "Join" button are two toggle controls: "Mic on" and "Camera on". At the bottom of the interface are two buttons: a grey "Back" button and a red "Remove from Favorites" button. Numbered callouts (1-6) in yellow circles with dashed arrows point to the following elements: 1 points to the "Mic on" toggle, 2 points to the "Camera on" toggle, 3 points to the "Join" button, 4 points to the "Remove from Favorites" button, 5 points to the "Remove from Favorites" button, and 6 points to the "Back" button.

# Dial Out to H.323 or SIP Endpoints

1

Click to open the Dial Out page

2

Click to enter an H.323 or SIP address (phone number)

3

Click to call (Emergency calls are not supported)

VidyoConnect™

Home

Search

Contacts

Rooms

Dial out

Settings

Help

20

a	b	c	d	e	f	
g	h	i	j	k	l	
m	n	<b>O</b>	p	q	r	
s	t	u	v	w	x	
y	z	0	1	2	3	
4	5	6	7	8	9	
Clear all	ABC	-	-	@	.	⌫


Mic on

Camera on

Call

Emergency Calls are not supported

# Meeting Controls (In-Call Experience)



## View Participants In-Call

1


Click to view the conference participants

2

View the display names of the participants


3

Click to pin a participant to a specific screen regardless of loudest speaker



# Control Camera In-Call

- 1 Click to access local camera controls
- 2 Use arrows to navigate up, down, left, and right
- 3 Use zoom buttons to zoom in and out
- 4 Click to return to main UI page





1

Click to access important contact names, numbers, emails, and web pages

2

View help details

3

Admins can customize text and upload a new help.html file in the Admin UI

The screenshot shows the VidyoConnect user interface. On the left is a vertical navigation menu with the following items: VidyoConnect logo, Home (house icon), Search (magnifying glass icon), Contacts (two people icon), Rooms (document icon), Dial out (grid of dots icon), Settings (gear icon), and Help (question mark icon). A yellow circle with the number '1' is positioned over the Help icon, with a dotted arrow pointing to it. A second yellow circle with the number '2' is positioned over the main content area, with a dotted arrow pointing to the top of the 'Healthcare team' page. The page content includes a large heading 'Healthcare team' and three paragraphs of text, each starting with 'Healthcare is a team effort. Each healthcare provider is like a member of the team with a special role. Some team members are doctors or technicians who help diagnose disease. Others are experts who treat disease or care for patients' physical and emotional needs.'

2

## Healthcare team

Healthcare is a team effort. Each healthcare provider is like a member of the team with a special role. Some team members are doctors or technicians who help diagnose disease. Others are experts who treat disease or care for patients' physical and emotional needs.

Healthcare is a team effort. Each healthcare provider is like a member of the team with a special role. Some team members are doctors or technicians who help diagnose disease. Others are experts who treat disease or care for patients' physical and emotional needs.

Healthcare is a team effort. Each healthcare provider is like a member of the team with a special role. Some team members are doctors or technicians who help diagnose disease.

Others are experts who treat disease or care for patients' physical and emotional needs.

1